MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICE RECEIVED

By Carol Day at 8:59 am, Mar 10, 2015

STATE PUBLIC HEALTH LABORATORY BREATH ALCOHOL PROGRAM

INTOX EC/IR II MAINTENANCE REPORT

REPORT #3

INTOX EC/IR II	MAINTENANCE AS	reventive mainte	mance check (not to	exceed 35
Complete this report at the time of days). Complete this report whenever	the regular monthly b	erviced or xepair	red and whenever it	is placed
days). Complete this report whenever into service. Retain the original a	and send a conv within	15 days to the f	Breath Alcohol Progr	cam, Dass
into service. Retain the original	NAME OF AGENCY		DATE OF INSPECTION	
INTOX EC/IR II SH	MARYLAND HEIGHTS POLICE		12/31/2014	
12825			TIME OF INSPECTION	
LOCATION OF INSTRUMENT (STREET AND CITY	, SHTS	,	00:13 CST	
11911 DORSETT RD. MARYLAND HEIC CHECKLIST: Place a mark in the box	hy each item if found	to be satisfacto	ory or is operating	.within
CHECKLIST: Place a mark in the box established limits. (Write in obse	rved values where deter	rmined). Unmark	ed items must be co	ITecred
established limits. (Write in obse	**************************************			
before using instrument.				
X DIAGNOSTIC RECORD	n n	X COZ CHECK		Annual Comments
X BLANK CHECK		X FLOW CHECK		
X FC 1 TEMP	X FCB CHECK			
X SRC TEMP	. Z		ЭK	
X SRC TEMP X CRC COMP CHECK X DET TEMP X CRC CAL CHECK				
X BT TEMP	100			
X STD 2 TEMP	PRINT TEST			
X STD 2 YEAR				
BREATH ANALYZER ACCURACY STAND	ARDS			מ
BREATH ANALIZER ACCORDED STATE		·	THANOL-GAS MIXTUR	B
X SIMULATOR SOLUTION	LABORATOR1ES L	OT# 14030		DATE 01/20/2016
I A IS WIADWO DOLLDIA	TO YOU THE THE THE		SIMULATOR EXP D	ATE
X SIMULATOR TEMP (34°C ±0.2°C	· 1		05/13/2015	
TY CALIBRATION CHECK - (ONLY OF	NE STANDARD IS TO BE	HED PER MAIN	THE PARTY AND COME OF SELECTION	of the standard value
Run three tests using a star	ndard solution. All	three tests m	ust be within +5%	ndard solution being
I and must have a spicac of .	005 or less. Mark t	the box corresp	ontorna to one ace	
used. (PRINTOUT ATTACHED)			1 7 170	
PERD RETURN 0.095% AND 0.105% INCLUSIVE				
X 0.10% STANDARD - MUST READ BETWEEN 0.076% AND 0.084% INCLUSIVE 0.08% STANDARD - MUST READ BETWEEN 0.038% AND 0.042% INCLUSIVE 0.04% STANDARD - MUST READ BETWEEN 0.038% AND 0.042% INCLUSIVE				
0.04% STANDARD - MUST REA	D BETWEEN 0.038% ANI	D 0.042% INCLUS	⊃ ⊤ Α 1 2	
			TEST 3 0.096	3 g/210L
TEST i * 0.095 g/210L	TEST 2 19 0.096	3/510P	TAN MUS VACIS MATERIA	TENANCE REPORT:
TEST 1 * 0.095 g/210L INDICATE THE NUMBER OF BREATH	TESTS IN THE FOLLOW	WING RANGES SIN	NCE THE DARK MAIN.	
INDICATE THE ROLL OF			.1519 5	OVER .19 0
REFUSALS 2 004 0	.0509 2	1014 3	TO . LD THOMPHUDAT	
REFUSALS 2 004 0 LIST ANY NEW PARTS AND DESCRIBE ANY A	LITERATION OR MODIFICATION	THAT WAS MADE TO	KERIOUP IND THRIVOUGHT	-
1 1 1 1 1 1 5		E december	,	
I MELIC NUMBER		E GEORGIANA.		
MEE-75 DHSS SI		E MECHANIAN.		
MEE-15 DH 53 7		. ABERRALATI	· .	
NEE-15 DH 53 7		, ABOURNALLY		
THEPECTING OFFICER		PRINT FULL NAME		
THEPECTING OFFICER	des Dof	PRINT FULL NAME SCHNURR, JAN	MES	
THE PECTING OF TEER		PRINT FULL NAME	MES	

Breath Alcohol Program, Missouri Department of Health and Senior Services,

Southeast District Office, 2875 James Blvd, Poplar Bluff, MO 63901

05/19/2016

(314)298-8700

240250

CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 14030 of
Alcohol Reference Solution for Simulator were analyzed by
gas chromatography on January 22, 2014, using a Perkin Elmer Gas
Chromatograph Autosystem XL S/N: 610N9030209, and found to contain
0.1215% (w/vol) ethyl alcohol. The expiration date for this lot
number is January 20, 2016 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:
Testing was conducted using Carilliant Reference Standard lot number EN122211-02 whose values are traceable to NIST:
All balances are calibrated annually by an outside agency using NIST traceable weights.
Calibration verification is done prior to each use utilizing NIST traceable weights.

14-10 Tab-035 12M

EXPIRES 5/19/2016



STATE OF MISSOURI

DEPARTMENT OF HEALTH AND SENIOR SERVICES BREATH ALCOHOL PROGRAM

ERMIT IL BULL

AMES R SCHO

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

DATAMASTER, INTOX EC/IR II

for the determination of the alcoholic content of blood from a sample of expired eir. Permit issued under the previsions of sections 577.02b through 577.041, RSMo and 306.111 through 306.119 RSMo.

NUMBER 240250	DATE 5/19/2014
Gal Vaclady	DIRECTOR OF STATE PUBLIC HEALTH LABORATORY

DIRECTOR OF DEPARTMENT OF HEALTH AND SENIOR SERVICES

764年5