

MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY BREATH ALCOHOL PROGRAM

DATAMASTER MAINTENANCE REPORT

REPORT #6

			1 / 11 105 doug)			
Complete this report at	the time of the regular monthly penever the instrument is service	oreventive maintenance	check (not to exceed 35 days). ever it is placed into service.			
Retain the original and	send a copy within 15 days to the	e Breath Alcohol Progra	am, DHSS.			
DATAMASTER SN	NAME OF ACENCY		DATE OF INSPECTION OZ/ZC/ZO/ZO/4 TIME OF INSPECTION			
940258 LOCATION OF INSTRUMENT (ST	BEET AND CITY)	AIC OF MILIO	I TIME OF INSPECTION			
Sag /a	A. A. COUTS TA	cy	0628 HOURS			
CHECKLIST: Place a m	lark in the box by each item it fou	nd to be satisfactory or i	f operating within established limits. (Write in observed values			
where determined.) Un	marked items must be corrected	before using instrumen	t.			
DIAGNOSTIC CHI	ECK (PRINTOUT ATTACHED)	DATE A	ND TIME (from printout) 02/20/14 0628 Hours			
COMPUTER		DET	ECTOR			
PROGRAM		FILT	<u>ERS</u>			
HEATERS SAN	MPLE CHAMBER4	9 °C ☑QUA	RTZ STANDARD			
FLOW DETEC	TOR	CAL	IBRATION			
PUMP HIGH S	PEED	PRIN	NTER			
M INDICATOR LIGH	TS					
SIMULATOR SOL	UTION SUPPLIER _ Gut	Н	LOT # 130 10 EXP. DATE 01 09 2015			
SIMULATOR TEM	IP (34°C ± 0.2°C)34	. O °C SIMULATOR	R SN <u>SD2505</u> EXP. DATE <u>01/23/2015</u>			
CALIBRATION CH	ECK - (ONLY ONE STANDARI	IS TO BE USED PER	MAINTENANCE REPORT)			
less. Mark the box	Run three tests using a standard solution. All three tests must be within ±5% of the standard value and must have a spread of .005 or less. Mark the box corresponding to the standard solution being used. (PRINTOUT ATTACHED)					
	DARD - MUST READ BETWEEN		· ·			
0.080% STANI	DARD - MUST READ BETWEEN	N 0.076% AND 0.084% I	INCLUSIVE			
☐ 0.040% STANI	DARD - MUST READ BETWEEN		INCLUSIVE			
TEST 1 . O	99 % TEST 2	100 %	TEST 3 . 100 %			
	TEST (PRINTOUT ATTACHED)	9				
INDICATE THE NUM	BER OF BREATH TESTS IN TH SELF-ADMINISTERED TESTS)	IE FOLLOWING RANG	ES SINCE THE LAST MAINTENANCE REPORT:			
	004) (.0509)	(.1014)	(.1519) OVER .19			
LIST ANY NEW PARTS AND DE	SCRIBE ANY ALTERATION OR MODIFICATION ARY).	N THAT WAS MADE TO RESTORE	THE INSTRUMENT TO OPERATE SATISFACTORILY AND WITHIN ESTABLISHED LIMIT			
		a a				
INSPECTING OFFIC	ER		PRINT EULL NAME			
SIGNATURE	Colon de		T SERNIE CLAYTON			
TYPE II PERMIT NUMBER/EXI	12 27 2014	A.	TELEPHONE NUMBER (636) 528-4725			
RETURN COMPLETED RI	PORT TO THE: Breath Alco	phol Program, MO Depa	rtment of Health and Senior Services, Southeast District Office			
	2875 Jame Poplar Bluf	s Blvd. f, MO 63901				



CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 13010 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on January 14, 2013, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1218% (w/vol) ethyl alcohol. The expiration date for this lot number is January 9, 2015 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN122211-02 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.

Face This Side Down - This Edge In First

BAC DataMaster **Evidence Ticket**

BAC DATAMASTER SERIAL NUMBER 948258

RREST TIME: 86:98

SUBECT NAME:

STATE/D.L.: MO/123456789 MDE SIANAL

ARREATING OFFICER

CLRTON/BERNIE FILER I.D.: 1 ETING OFFICER: 31

CLAYTON/BERNIE OFFICER I.D.: 131

PERMIT NUMBER: 229414

EXPIRATION DRITE: 12/27/14 MISCELLAMEDUS DATA: BILLA 35 INTERNAL STANDARD
- SUBJECT SAMPLE
BLANK TEST 医牙术 古

8888 8888

BAC DataMaster Evidence Ticket

Face This Side Down - This Edge In First

DAC DATAMASTER SERIAL NUMBER 940258 82/28/14

HREST IME: 66:88 SECT MAR

RFICTEST

DOB: GI/GI/GI

STATE/D.L.: MO/123456789 ARRESTING OFFICER:

CLATOVBENIE OFFICER I.D.: 131 TEXTING OFFICER:

CHIONERNIE

PERMIT NUMBER: SERIA OFFICER L.D.: 131

EXPIRATION WATE: 12/27/14 - BRENTH ANALYSIS

RADIO INTERFERENCE INTERNAL STANDARD RINK TEST

Operator Signature

Operator Signature_

Face This Side Down - This Edge In First

BAC DataMaster Evidence Ticket

STATE OF MISSOURI BAC DATAMASTER SERIAL MUMBER 940258 02/20/14

ESTIME OFFICER:

OFFICER I.D.: 131 PERMIT NUMBER: 228414

EXPIRATION DATE: 12/27/1 MISCELLAMEOUS DATA:

--- SUPERVISOR MORE

13								29.5	
		S	8	8	S	Š	8	S	
	S S	VERITIE	er.			995		99.	
			発見を		EXTERNE STRUKE		1	BLRK ROT	

Face This Side Down - This Edge In First

BAC DataMaster **Evidence Ticket**

ENC DATHMASTER SERIAL NUMBER 940258 STATE OF MISSOURI 3 3

名誉 COMPUTER:

FS PROGRAM (64-87-2889):

HENTERS SAMPLE CHRMBER:

ij GK St

FLOW DETECTOR:

ë Se

E, A B

DETECTOR:

K LIES.

WHITE STANDED

CHLERRILOM

ă

844% (0*4*-1.70123455789:;<=>?*RDCDEFG HIJKLMYOPORSTUWWXYZIJ~Jabcdefahijklmmo **米加尼** 四 Parsturumxsz(|)+4

Operator Signature

Operator Signature_

State of Missouri DEPARTMENT OF HEALTH



PERMIT TYPE II



BERNIE CLAYTON

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field repairs, and operate the following breath analyzer(s):

DATAMASTER

for the determination of the alcoholic content of blood from a sample of expired (alveolar) air. Issued under the provisions of sections 577.020 through 577.041, RSMo 1986.

12/27/2012	Wante
Number 220414	Director of State Public Health Laboratory
Expires 12/27/2014	Lail Vacherly ACTING DIRECTOR
MO 580-0771 (7-88)	Director, Department of Health Lab. 4 (R7-88)