

## MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY

BREATH ALCOHOL PROGRAM

#### DATAMASTER MAINTENANCE REPORT

RECEIVED . By Carol Day at 9:32 am, Aug 06, 2013

REPORT#6

·								
Complete this report at the time of the regular monthly preventive maintenance check (not to exceed 35 days).  Complete this report whenever the instrument is serviced or repaired and whenever it is placed into service.  Retain the original and send a copy within 15 days to the Breath Alcohol Program, DHSS.								
DATAMASTER SH NAME OF AGENCY					DATE OF INSPECTION	1		
2012		St. Joseph Police	ce Department		08/01/2013			
l .	том о <i>г</i> інстримен Faraon St. Jose	T (STREET AND CITY)			TIME OF INSPECTION 2013			
			ch item if found to be s	satisfactory or if operal	· · · · · · · · · · · · · · · · · · ·	nits. (Write in observed values		
CHECKLIST: Place a mark in the box by each item if found to be satisfactory or if operating within established limits. (Write in observed values where determined.) Unmarked items must be corrected before using instrument.								
Ø	DIAGNOSTIC CHECK (PRINTOUT ATTACHED)			DATE AND TIME	DATE AND TIME (from printout) 08/01/2013 2020			
	COMPUTE	R		DETECTOR				
	N PROGRAM		,	FILTERS				
	HEATERS	SAMPLE CHAMBER 48	oc	QUARTZST	ANDARD			
	FLOW DET	ECTOR	as more		N			
	PUMP HIGI	H SPEED						
×	INDICATOR LIC	GHTS						
$\boxtimes$	SIMULATOR S	OLUTION SUPPLIER R	epco Marketing	LOT#	12001 EXF	RDATE 04/02/2014		
$\boxtimes$	SIMULATOR T	EMP (34'C ± 0.2'C) 34.	0 ос	SIMULATOR SN SD	2278 EXF	P, DATE 12/17/2013		
$\nabla$	CALIBRATION (	CHECK - (ONLY ONE S	TANDARD IS TO BE	USED PER MAINTE	NANCE REPORT)			
						and bear a surred of OOE on		
	Run three tests using a standard solution. All three tests must be within ±5% of the standard value and must have a spread of .005 or less. Mark the box corresponding to the standard solution being used. (PRINTOUT ATTACHED)							
, 	KT3							
ř	O.100% STANDARD - MUST READ BETWEEN 0.095% AND 0.105% INCLUSIVE 0.080% STANDARD - MUST READ BETWEEN 0.076% AND 0.084% INCLUSIVE							
	0.040% STANDARD - MUST READ BETWEEN 0.076% AND 0.044% INCLUSIVE							
TES	T 1098		TEST 2100		TEST 3 .099			
		I. TEST (PRINTOUT AT						
_			· · · · · · · · · · · · · · · · · · ·	UINO DANOES SINO	E THE LAST MAINTE	ANCE DEBORT		
		SELF-ADMINISTERED		VIIVO KANGES SINC	C INCLAST WAINTON	ANGE REFORT.		
				1		1		
REF	USALS 0	(004) <sub>0</sub>	(.0509)	(.1014) 0	(.1509) 0	(OVER .19) <sub>1</sub>		
			MODIFICATION THAT WAS MA	DE TO RESTORE THE INSTRU	MENT TO OPERATE SATISFACTO	RILY AND WITHIN ESTABLISHED LIMITS		
(USE OTHER SIDE IF NECESSARY).								
	`							
NEDECTING OFFICED								
INSPECTING OFFICER SIGNATURE PRINT FULL NAME								
1	1/1//		5//	Scott Gar	Scott Gary Brad Kerns			
	PEMIT NUMBERIEX		2711	TELEPHONE	TELEPHONE NUMBER			
2017	79 08/03/2013	220427 12/27/2014			816-271-5359			
RETURN COMPLETED REPORT TO THE: Breath Alcohol Program, MO Department of Health and Senior Services, Southeast District Office 2875 James Blvd.								

Poplar Bluff, MO 63901

#### CERTIFICATE OF ANALYSIS

MANUFACTURER AND SUPPLIER: RepCo Marketing, Inc.
LOT NUMBER: 12002
EXPIRATION DATE: August 29, 2014 at 11:59 p.m.

RepCo Marketing, Inc. certifies the following:

RepCo Marketing, Inc. manufactured, tested and supplied Lot Number 12002 of Alcohol Certified Solution for simulators. Random samples of said lot number were analyzed by an independent laboratory utilizing a gas chromatograph gms/dl +/-.003 gms/dl wt/vol ethanol (95% and found to contain .1209 Confidence).

The alcohol and distilled water used in the solution were found to be free of any interferring substance.

This solution will produce a vapor alcohol value of \_.100 +/-3% gms/210L Breath when heated to 34 Degrees Celsius +/-0.2 Degrees Celsius in a simulator (95% Confidence).

The date of manufacture for this lot number is August 30, 2012 The expiration date for this lot number is \_\_\_\_Awast 29, 2014 11:59 p.m.

This document is a true representation of the original Certificate of Analysis.

Cecil B. Garner, President

RepCo Marketing, Inc.

Form RM 02

## FACE THIS SIDE DOWN - THIS EDGE IN FIRS

### **BAC DataMaster**

**Evidence Ticket** 

STATE OF MISSOURI ST. JOSEPH POLICE DEPARTMENT

BAC DATAMASTER SERIAL NUMBER 201206 08/01/13 20:20

--- DIAGNOSTIC CHECK ---

COMPUTER:

UKAY

PROGRAM (04-07-2009):

OKAY

HEATERS

SAMPLE CHAMBER:

48c

FLOW DETECTOR:

OKAY

PUMP

HIGH SPEED:

OKAY

DETECTOR:

OKAY

FILTERS:

OKAY

QUARTZ STANDARD:

OKAY

CALIBRATION:

OKAY

PRINTER TEST

!"#\$%&'<)\*+.-./0123456789#;<\\?@ABCDEFG HIJKLMNOP@RSTUVWXYZ[\]^\_\abcdefghijklmno pqrstuvwxyz{|}}~

#### **OPERATOR SIGNATURE**

Card Stock No.

60021

REORDER ALL SUPPLIES FROM N.P.A.S. P.O. BOX 1435, MANSFIELD, OH 44901

#### FACE THIS SIDE DOWN - THIS EDGE IN FIRST

### **BAC DataMaster**

**Evidence Ticket** 

STATE OF MISSOURI ST.JOSEPH POLICE DEPARTMENT

BAC DATAMASTER SERIAL NUMBER 201206 08/01/13

TESTING OFFICER:
GARY/SCOTT
OFFICER I.D.: 1013
PERMIT NUMBER: 220179
EXPIRATION DATE: 08/03/14
MISCELLANEOUS DATA:

#### --- SUPERVISOR MODE ---

BLANK TEST	.000	2 <b>0:</b> 25
INTERNAL STANDARD	VERIFIED	20:25
EXTERNAL STANDARD	.098	20:25
BLANK TEST	`.000	20:26
EXTERNAL STANDARD	.100	20:26
BLANK TEST	. 000	20:27
EXTERNAL STANDARD	.099	20:27
BLANK TEST	.000	20:28

√N = 3 | SIM. = .1 | AVG. = .099

**OPERATOR SIGNATURE** 

Card Stock No.

REORDEH ALL SUPPLIES FROM N.P.A.S. RO. BOX 1435, MANSFIELD, OH 44901

### FACE THIS SIDE DOWN - THIS EDGE IN FIRST

#### **BAC DataMaster**

**Evidence Ticket** 

STATE OF MISSOURI ST.JOSEPH POLICE DEPARTMENT

BAC DATAMASTER SERIAL NUMBER 201206 08/01/13

ARREST TIME: 19:00 SUBJECT NAME: K DOB: 09/09/09 SEX: M STATE/D.L.: MO/9 ARRESTING OFFICER: ΚJ OFFICER I.D.: 8J TESTING OFFICER: GARY/SCOTT OFFICER 1.D.: 1013 PERMIT NUMBER: 220179 EXPIRATION DATE: 08/03/14 MISCELLANEOUS DATA: RFI TEST

#### --- BREATH ANALYSIS ---

<sup>1</sup>BLANK TEST INTERNAL STANDARD RADIO INTERFERENCE .000 20:13 VERIFIED 20:13

**OPERATOR SIGNATURE** 

Card Stock No. 60021

REORDER ALL SUPPLIES FROM N.P.A.S. P.O. BOX 1435, MANSFIELD, OH 44901

# State of Missouri DEPARTMENT OF HEALTH



# PERMIT TYPE II



#### **BRAD M KERNS**

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field repairs, and operate the following breath analyzer(s):

#### DATAMASTER; ALCO-SENSOR IV W/PRINTER

for the determination of the alcoholic content of blood from a sample of expired (alveolar) air. Issued under the provisions of sections 577.020 through 577.041, RSMo 1986.

12/27/2012	Wenter		
	Director of State Public Health Laboratory		
Number 220427	Dal Vooledy Acting Director		
Expires 12/27/2014	ACTING DIRECTOR		
EVS/1109	Director, Department of Health		
MO 580-0771 (7-88)	Lab. 4 (A7-88)		