

#### MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY BREATH ALCOHOL PROGRAM DATAMASTER MAINTENANCE REPORT

RECEIVED

By Carol Day at 4:03 pm, Sep 06, 2013

REPORT #6

Complete this report at the time of the region Complete this report whenever the instrum Retain the original and send a copy within	ient is serviced or repaired	l and whenever it is i	nlaced into convice	),	
NAME OF AGENCY 201237 Lee's Summit Police Department			DATE OF INSPECTION 09/05/2013		
LOCATION OF INSTRUMENT (STREET AND CITY)  10 NE Tudor Rd, Lee's Summit, MO 6	`	TIME OF INSPECTION 3:50 pm			
CHECKLIST: Place a mark in the box by ea where determined.) Unmarked items must	ich item if found to be satis	factory or if operating instrument.		nits. (Write in observed values	
DIAGNOSTIC CHECK (PRINTOUT AT	TACHED)	DATE AND TIME (	from printout) 09/05	/2013 1550 HRS	
☑ COMPUTER		☑ DETECTOR			
☑ PROGRAM	·. ·	☑ FILTERS			
☐ HEATERS SAMPLE CHAMBER _	<u>49.0</u> °C	Quartz Stan	IDARD		
☑ FLOW DETECTOR		☐ CALIBRATION			
PUMP HIGH SPEED		PRINTER			
INDICATOR LIGHTS					
SIMULATOR SOLUTION SUPPLIER	Guth Laboratories	LOT # _12	2040 EXF	P, DATE 03/07/2014	
SIMULATOR TEMP (34°C ± 0.2°C)	34.0 °C SIN	MULATOR SN	SD3148 EX	P. DATE 01/07/2014	
A CALIBRATION CHECK - (ONLY ONE S	STANDARD IS TO BE US	ED PER MAINTENA	ANCE REPORT)	****	
Run three tests using a standard solution. All three tests must be within ±5% of the standard value and must have a spread of .005 or less. Mark the box corresponding to the standard solution being used. (PRINTOUT ATTACHED)					
0.100% STANDARD - MUST READ 0.080% STANDARD - MUST READ 0.040% STANDARD - MUST READ	BETWEEN 0.095% AND ( BETWEEN 0.076% AND (	0.105% INCLUSIVE 0.084% INCLUSIVE	· ·	-	
TEST 1100	TEST 2 ₩ .100		TEST 3 🕶 .098	****	
PERFORM R.F.I. TEST (PRINTOUT AT	rached)				
INDICATE THE NUMBER OF BREATH TESTS IN THE FOLLOWING RANGES SINCE THE LAST MAINTENANCE REPORT: (DO NOT INCLUDE SELF-ADMINISTERED TESTS)					
REFUSALS 0 (004) 0	(.0509) 0 (.	1014) 1	(.1519) 0	OVER .19 1	
IST ANY NEW PARTS AND DESCRIBE ANY ALTERATION OR USE OTHER SIDE IF NECESSARY).  None	MODIFICATION THAT WAS MADE TO	RESTORE THE INSTRUMEN	T TO OPERATE SATISFACTOR	RILY AND WITHIN ESTABLISHED LIMITS	
				·	
NSPECTING OFFICER	Separate extension	BASA BERTHAN			
Chandes Kudinales		PRINT FULL NAME MPO Jennifer L Kudzinski			
PPE II JERMIT NUMBER/EXPIRATION DATE 06/12/2015			TELEPHONE NUMBER (816) 969-4150		
Breath Alcohol Program, MO Department of Health and Senior Services, Southeast District Office 2875 James Blvd. Poplar Bluff, MO 63901					

# FACE THIS SIDE DOWN - THIS EDGE IN FIRST

# **BAC DataMaster**

**Evidence Ticket** 

STATE OF MISSOURI LEES SUMMIT POLICE DEPARTMENT

BAC DATAMASTER SERIAL NUMBER 201237 09/05/13

TESTING OFFICER: KUDZINSKI/JENNIFER OFFICER 1.D.: 0524 PERMIT NUMBER: 230121 EXPIRATION DATE: 06/12/15 MISCELLANEOUS DATA:

#### --- SUPERVISOR MODE ---

BLANK TEST	.000	15:38
INTERNAL STANDARD	VERIFIED '	15:38
EXTERNAL STANDARD	. 100	15:39
BLANK TEST	.000	15:40
EXTERNAL STANDARD	. 100	15:40
BLANK TEST	.000	15:41
EXTERNAL STANDARD	.098	15:41
BLANK TEST	.000	15:42

N = 3 SIM. = .1 AVG. = .0993

OPERATOR SIGNATURE Card Stock No. 60021

REORDER ALL SUPPLIES FROM N.P.A.S. P.O. BOX 1435, MANSFIELD, OH 44901

## FACE THIS SIDE DOWN - THIS EDGE IN FIRST

### **BAC DataMaster**

**Evidence Ticket** 

STATE OF MISSOURI LEES SUMMIT POLICE DEPARTMENT

BAC DATAMASTER SERIAL NUMBER 201237 09/05/13 15:50

--- DIAGNOSTIC CHECK ---

COMPUTER:

ÖKAY

PROGRAM (04-07-2009):

OKRY

HEATERS

SAMPLE CHAMBER:

49c

FLOW DETECTOR:

OKAY

PUMP

HIGH SPEED:

OKAY

DETECTOR:

OKAY

FILTERS:

OKAY

QUARTZ STANDARD:

CALIBRATION:

OKAY

ÖKAY

#### PRINTER TEST

!"#\$%&'()\*+,~./0123456789:;<=>?@ABCDEFG HIJKLMNOPQRSTUVWXYZI\l^\_'abcdefghijklmno pqrstuvwxyz{|}\*~

OPERATOR SIGNATURE

REORDER AL SUPPLIES FROM N.P.A.S. P.O. BOX 1435, MANSFIELD, OH 44901

FACE INIO SIVE DOWN - INIO EDUC IN FIRST

### **BAC DataMaster**

**Evidence Ticket** 

STATE OF MISSOURI LEES SUMMIT POLICE DEPARTMENT

BAC DATAMASTER SERIAL NUMBER 201237 09/05/13

ARREST TIME: 15:15 SUBJECT NAME: TEST DOB: 10/10/90 STATE/D.L.: MO/11111111 ARRESTING OFFICER: KUDZINSKI/JENNIFER OFFICER I.D.: 0524 TESTING OFFICER: KUDZINSKI/JENNIFER OFFICER I.D.: 0524 PERMIT NUMBER: 230121 EXPIRATION DATE: 06/12/15

--- BREATH ANALYSIS ---

BLANK TEST

.000

15:48

INTERNAL STANDARD RADIO INTERFERENCE

MISCELLANEOUS DATA:

VERIFIED

15:48

OPERATOR SIGNATURE

Card Stock No.

REORDER ALL SUPPLIES FROM N.P.A.S. P.O. BOX 1435, MANSFIELD, OH 44901



#### STATE OF MISSOURI

DEPARTMENT OF HEALTH AND SENIOR SERVICES BREATH ALCOHOL PROGRAM



# JENNIFER KUDZINSKI

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

## DATAMASTER, INTOXILYZER 5000, ALCO-SENSOR IV W/PRINTER

	rough 577.041, HSMo and 306.1		SMo.
DATE	06/12/2013	, 	Wenter
_,	230121		DIRECTOR OF STATE PUBLIC HEALTH LABORATORY  LIGHT VOOLENCY
HOMBETT -	06/12/2015	• . •	Acting Director
EXPIRES =		_	DIRECTOR OF DEPARTMENT OF HEALTH AND SENIOR SERVICES

MC 680-0771 (6-10)

Reference Standard lot number FNI22211-02 whose an outside agency using NIST traceable weights.

GUTH LABORATORIES,

Pauley, President

test interfering substances.

alcohol and water used in this solution were analysis instrument reading of 0.100 g/210L +/- 3%

.2°C,

this solution will give a

breath alcoho

operating

calibrated Simulator,

0.1211% (w/yol) ethyl alcohol.

Alcohol Reference

Solution

for Simulator were

Sursn

Perkin analyzed

Elmer ₽

found

contain

chromatography on March

XI S/N: 610N9030209, and

ģ



Alcohol Reference

Solution

Simulator