

MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY BREATH ALCOHOL PROGRAM DATAMASTER MAINTENANCE REPORT

RECEIVED

By Carol Day at 10:05 am, Feb 03, 2014

REPORT #6

Complete this report at the time of the regular Complete this report whenever the instrume Retain the original and send a copy within 1.	nt is serviced or repaire	ed and whenever it is	placed into serv	days). rice.			
DATAMASTER SN NAME OF AGENCY Chillicothe Po				DATE OF INSPECTION 01/28/2014			
LOCATION OF INSTRUMENT (STREET AND CITY) 613 Walnut Street, Chillicothe, Missouri 64601			TIME OF INSPECTION 0:19 am				
CHECKLIST: Place a mark in the box by each where determined.) Unmarked items must be			g within establisl	ned limits.	(Write in observ	ved value:	
☑ DIAGNOSTIC CHECK (PRINTOUT ATT		E AND TIME (from printout) 01/28/2014 00:34					
☑ COMPUTER	✓ DETECTOR	ETECTOR					
✓ PROGRAM	✓ FILTERS	LTERS					
✓ HEATERS SAMPLE CHAMBER	✓ HEATERS SAMPLE CHAMBER						
✓ FLOW DETECTOR ✓ CALIBRATION							
✓ PUMP HIGH SPEED ✓ PRINT			ER				
✓ INDICATOR LIGHTS							
SIMULATOR SOLUTION SUPPLIER Guth Laboratories LOT # 13210 EXP. DATE 07/29/2015						15	
SIMULATOR TEMP (34°C ± 0.2°C)	34.0 °C S	IMULATOR SN	SD2270	_ EXP. DA	TE 12/06/20	14	
CALIBRATION CHECK - (ONLY ONE ST	ANDARD IS TO BE U	SED PER MAINTEN	ANCE REPORT				
CALIBRATION CHECK - (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT) Run three tests using a standard solution. All three tests must be within ±5% of the standard value and must have a spread of .005 or							
less. Mark the box corresponding to the standard solution being used. (PRINTOUT ATTACHED)							
0.100% STANDARD - MUST READ BETWEEN 0.095% AND 0.105% INCLUSIVE 0.080% STANDARD - MUST READ BETWEEN 0.076% AND 0.084% INCLUSIVE							
0.040% STANDARD - MUST READ B							
TEST 1 ▼ .096		TEST 3 • .098					
PERFORM R.F.I. TEST (PRINTOUT ATTA	ACHED)						
INDICATE THE NUMBER OF BREATH TEST (DO NOT INCLUDE SELF-ADMINISTERED	TS IN THE FOLLOWIN TESTS)	IG RANGES SINCE 1	THE LAST MAI	NTENANO	E REPORT:		
REFUSALS 0 (004) 6 ((.0509)	(.1014) 0	(.1519)	1	OVER .19	0	
LIST ANY NEW PARTS AND DESCRIBE ANY ALTERATION OR MO (USE OTHER SIDE IF NECESSARY).	DDIFICATION THAT WAS MADE 1	TO RESTORE THE INSTRUMEN	IT TO OPERATE SATIS	FACTORILY A	 ND WITHIN ESTABLIS	SHED LIMITS	
Guth Laboratories Inc., .100 Solution, Lo	t # 13210 , Bottle # 4	86. Expires 07/29/2	015 at 2359 h	ours			
		, , , , , , , , , , , , , , , , , , , ,					
INSPECTING OFFICER SIGNATURE		PRINT FULL NAM					
· Tokent I Helder	to resolve the resolvent to the second secon	Robert L. Heldenbrand					
TYPE II PERMIT NUMBER/EXPIRATION DATE 220161 07/06/20		TELEPHONE NUMBER (660) 646-2121					
RETURN COMPLETED REPORT TO THE: Breath Alcohol Program, MO Department of Health and Senior Services, Southeast District Office							

Poplar Bluff, MO 63901



CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 13210 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on July 31, 2013, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1216% (w/vol) ethyl alcohol. The expiration date for this lot number is July 29, 2015 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN122211-02 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.

Face This Side Down – This Edge In First

BAC DataMaster Evidence Ticket

STATE OF MISSOUR! CHILLICOTHE POLICE DEPRETAENT

BAC DATAPASTER SERIAL NUMBER 958415 01/28/14

ESTRE OFFICER

HELDETHRAND/RAL

OFFICER 1.D. # 11
PERMIT NUMBER: 228161
PATENTON DRIE: 37/26/14
MISCELLAYGOUS DRIE: 37/26/14

8 8 8 8 8 8 8 8 8 8 8 8 8 8 8 8 8 8 8	
	0 00 00 5 05 60 0 00 65
HATERAL STRUCKS EMTERNAL STRUCKS BLANK LST EMTERNAL STRUCKS BLANK LEST MANNEN BLANK LEST MANNEN BLANK LEST	EXTENDED STRUGGED BLANK TEST

(O H 22

Face This Side Down – This Edge In First

BAC DataMaster Evidence Ticket

STATE OF MISSOURI CHILLICOTHE POLICE DEPRETMENT

BAC DATAMANTER SERIAL NUMBER 958415 81/28/14

J MARGARA MARGARA

MENTES STATE CHARES

U Gh Th

FLOW METECTURY

MILCIPA

QUARTZ STRUMRDI

!"#\$X&'()**+../B123456789::<->?eABCDEF6 HIJKLMNOPGRSTUVWXYZI\J\] abodefshiklmo 李明"生意思,是这种思想的是一个多个

Operator Signature

2208-02

2208-02

Operator Signature

Face This Side Down - This Edge In First

BAC DataMasterEvidence Ticket

STATE OF MISSOURI
CHILLICOTHE POLICE DEPARTMENT

BAC DATAMASTER SERIAL NUMBER 950415 01/28/14

ARREST TIME: 80:01
SUBJECT NAME:
RFI/TEST
DOB: 01/01/01 SEX: M
STATE/D.L.: MO/123456789
ARRESTING OFFICER:
HELDENBRAND
OFFICER I.D.: 11
TESTING OFFICER:
SAME
OFFICER I.D.: 11
PERMIT NUMBER: 220161
EXPIRATION DATE: 07/06/14
MISCELLANEOUS DATA:
JAN14 MRINTENANCE

--- BRESTH PARLYSIS ---

RADIO IMTERFERENCE

Operator Signature_

2208-02