

MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES

STATE PUBLIC HEALTH LABORATORY BREATH ALCOHOL PROGRAM

RECEIVED

By Carol Day at 11:54 am, Mar 30, 2015

INTOX EC/IR II MAINTENANCE REPORT

REPORT #3

	I MAINIENANCE			REPORT #1
Complete this report at the time of the regular monthly preventive maintenance check (not to exceed 35 days). Complete this report whenever the instrument is serviced or repaired and whenever it is placed				
into service. Retain the original		in 15 days to the		
INTOX EC/IR II SN	NAME OF AGENCY		DATE OF INSPECTION	N .
12858	Holts Summit PD		03/27/2015	
LOCATION OF INSTRUMENT (STREET AND CIT			TIME OF INSPECTION	N
245 South Summit Dr Holts Sum			15:27 CDT	
CHECKLIST: Place a mark in the bo				
established limits. (Write in obs	erved values where de	etermined). Unmar	ced items must be	corrected
before using instrument. X DIAGNOSTIC RECORD				
		TWIGO? OUTOW		
X BLANK CHECK	944 mm.	X CO2 CHECK		
X FC 1 TEMP		X FLOW CHECK		
X SRC TEMP		X FCB CHECK		
X DET TEMP		X CRC COMP CHEC	CK	
X BT TEMP		X CRC CAL CHECK	K	
X STD 2 TEMP		X PRINT TEST		
Х ЕТН СНЕСК				
	A D D C			
BREATH ANALYZER ACCURACY STANI	ARDS	ET COMPAGE DE	NILLIAN CLASSICAL	TD D
SIMULATOR SOLUTION			THANOL-GAS MIXTU	
		LOT# AG414702		DATE 05/27/2016
SIMULATOR TEMP (34°C ±0.2°C) SIMULA	TOR S/N	SIMULATOR EXP	DATE
X CALIBRATION CHECK - (ONLY OR	E STANDARD IS TO B	E USED PER MAINT	ENANCE REPORT)	
Run three tests using a star				% of the standard value
and must have a spread of .005 or less. Mark the box corresponding to the standard solution being used. (PRINTOUT ATTACHED)				
X10.10% STANDARD - MUST REAL	BETWEEN 0.095% AN	D 0.105% INCLUSI	VE	
0.08% STANDARD - MUST REAL				
0.04% STANDARD - MUST REAL	BETWEEN 0.038% AN	D 0.042% INCLUSI	VE	
TEST 1 0.098 g/210L	TEST 2 0.098	g/210L	TEST 3 0.09	8 g/210L
INDICATE THE NUMBER OF BREATH	TESTS IN THE FOLLO	WING RANGES SINC	E THE LAST MAIN	TENANCE REPORT:
REFUSALS 1 004 0	.0509 0	.1014 0	.15~.19 0	OVER .19 0
LIST ANY NEW PARTS AND DESCRIBE ANY AL	ERATION OR MODIFICATION	THAT WAS MADE TO RES	STORE THE INSTRUMENT	TO OPERATE
SATISFACTORILY AND WITHIN ESTABLISHED	IMITS (USE OTHER SIDE I	F NECESSARY).		
INSPECTING OFFICER				
SIGNATURE		PRINT FULL NAME		
> Bund Kil		REID, BRYAN		
	ATION DATE	TELEPHONE NUMBER		
240316 / 07/	30/2016	(573)896-8431		
RETURN COMPLETED REPORT TO THE:				
Breath Alcohol Program, Mis		of Health and	Senior Service	es,
-				,
Southeast District Office, 2875 James Blvd, Poplar Bluff, MO 63901				



Airgas USA LLC (LAB)

3500 Bernard Street

St. Louis, Mo. 63103

Ph: (314) 533-3100

Fax: (314) 533-7328

Certificate of Analysis

<u>Customer Name</u> Intoximeters, Inc. 2081 Craig Road St. Louis, Mo 63146 Test Date: 30-May-2014

Lot # AG414702

Exp. Date 27-May-2016 Cyl. Type

Component Ethanol Nitrogen Certified Concentration 0.100 ± 2% BrAC (272 ppm)

Balance

Certification Traceable to N.I.S.T. RGM Ethanol Standards:

Serial No.	Concentration	<u>Serial No.</u>	<u>Concentration</u>
EB0010581	391.8 ppm	EB0010603	392.5 ppm
EB0010570	259.8 ppm	EB0010559	258.9 ppm
EB0010285	209.0 ppm	EB0010595	208.9 ppm
EB0010561	103.7 ppm	EB0010562	104.9 ppm
EB0010681	52.22 ppm	EB0010579	52.94 ppm

Analytical Method:

NDIR

Analyst:

Pod Marcala

ISO 17025:2005 A2LA accredited. Certificate Number 2989.01

U.S. Department of Commerce National Institute of Standards and Technology Material Measurement Laboratory Chemical Sciences Division Gaithersburg, MD 20899-8393

REPORT OF ANALYSIS

November 15, 2012

Recertification of Ten (50-400) μmol/mol Ethanol in Nitrogen RGM Standards for Airgas Mid-America, St. Louis, MO

Submitted to:

Mr. Randy Renner Airgas Mid-America 3500 Bernard Street St. Louis, MO 63103

Job No.: 13010

P.O. No.: 4501370879

Airgas Mid-America (AMA) submitted ten reference gas material (RGM) standards (nominal 50 to 400 μmol/mol ethanol in nitrogen) to the National Institute of Standards and Technology (NIST) for recertification. These standards were analyzed in 2009 by gas chromatography / flame ionization detector (GC/FID). This analytical technique is difficult due to the tailing of the ethanol peak (especially at high concentration) in the gas chromatogram. To avoid this, NIST recently developed a Fourier transform infrared (FTIR) method to analyze ethanol, as described below. This new technique was used in a bilateral comparison of 120 μmol/mol ethanol (nitrogen balance) with the Dutch Metrology Institute (Von Schwinden Laboratories). The relative difference between the NIST and VSL-assigned concentrations was 0.12 % [1]. Consequently, these standards were analyzed by FTIR.

The current FTIR analyzed concentrations are within the expanded uncertainty limits of the 2009 certified values. However, since FTIR is a superior analytical technique over GC/FID (for ethanol), then the current analyzed values are a better reflection of the true ethanol concentrations of the samples. These samples are recertified at the current analyzed concentration as shown in table 1. The recertification of these RGM standards was in accordance with the Gas Metrology Group Quality Manual (QM-III-839.03), TP 839.03.11B.

The uncertainties in the certified concentrations are lower than those reported in the 2009 certificate and are expressed as an expanded uncertainty $U = ku_c$ with u_c estimated from the experimental standard deviations and the coverage factor k equal to 2. The true concentration is asserted to lie within the interval defined by the certified value $\pm U$ with a level of confidence of approximately 95 % [2].

Analytical Method

Ethanol Analysis: The ethanol (EtOH) component of each submitted sample was analyzed by Fourier Transform Infrared (FTIR) (Nicolet Model Nexus 670, NIST # 593134) equipped with a 10 meter (m) folding path, quartz gas cell (Specac, Model Cyclone 10C, NIST# 623477) with potassium bromide windows and a mercury-cadmium-telluride (MCT) detector.

Calibration

Two calibration curves were developed using the primary gas standards listed in table 2. One curve for the range 24 to 150 μ mol/mol ethanol; the other for the range 150 to 440 μ mol/mol. The lower curve was linear, the higher curve was quadratic.

Traceability

The instrument was calibrated using NIST primary gas standards which are prepared gravimetrically from pure gases, verified, and periodically compared internationally with other national metrology institutes' primary gas standards. This assures traceability to the SI for gas standards certified by NIST. See table 2 for a listing of the primary gas standards used.

Certification Period

These certified values are valid for a period of four (4) years from the date of the report. Consequently, the expiration date for these standards is 11/15/2016. This report serves as the NIST Certificate of Certification for these materials.

References:

- 1. Analysis of an Ethanol Gas Standard from VSL as a Bilateral Comparison: ROA # 639.03-11-045 [12/7/2010]
- 2. Guide to the Expression of Uncertainty, ISBN 92-67-10188-9, 1st Edition, ISO, Geneva, Switzerland, 1993.

Other References: Notebook: FTIR LG#2 [ACD # 3668]; pp. 1-6

Original RGM Standard Certification: ROA # 839.03-09-041[02/10/2009]

Prepared by:

Lyn Gameson

Research Chemist

Reviewed by:

Franklin R. Guenther, Ph.D.

Michael C. Killy (Jossyr Gamera)

Leader,

Gas Metrology Group

Michael E. Kelley Research Chemist

Carlos A. Gonzales, Ph.D.

Chief

Chemical Sciences Division

Table 1: Recertification of AMA Ethanol in Nitrogen RGM Standards. The stated uncertainties are $\overline{\text{expanded (k = 2)}}$.

Cylinder Number	2009 Analysis EtOH (µmol/mol) ^a	2012 Analysis EtOH (μmol/mol) ^b	Current Certified Conc. (µmol/mol)	Pressure (MPa)
EB0010579	52.4 ± 1.1	52.94 ± 0.70	52.94 ± 0.70	13.6
EB0010681	53.0 ± 1.1	52.22 ± 0.70	52.22 ± 0.70	13.6
EB0010562	104.9 ± 2.1	104.9 ± 1.3	104.9 ± 1.3	11.2
EB0010561	101.9 ± 2.1	103.7 ± 1.3	103.7 ± 1.3	13.6
EB0010595	209.2 ± 4.2	208.9 ± 2.5	208.9 ± 2.5	13.3
EB0010285	208.9 ± 4.2	209.0 ± 2.5	209.0 ± 2.5	13.6
EB0010559	258.3 ± 5.2	258.9 ± 3.2	258.9 ± 3.2	12.4
EB0010570	258.4 ± 5.2	259.8 ± 3.2	259.8 ± 3.2	13.6
EB0010603	390.9 ± 7.8	392.5 ± 5.0	392.5 ± 5.0	9.5
EB0010581	391.5 ± 7.8	391.8 ± 5.0	391.8 ± 5.0	9.3

^a Original Certification (GC/FID) ^b Current Analysis (FTIR)

Table 2: Ethanol in Nitrogen NIST Primary Gas Standards. The stated uncertainties are expanded (k=2).

Cylinder Number	EtOH (µmol/mol)
AAL20661	24.06 ± 0.05
AAL20255	38.78 ± 0.08
ALM024319	54.95 ± 0.11
ALM040277	70.54 ± 0.14
ALM040295	100.7 ± 0.2
ALM009006	121.8 ± 0.2
ALM040288	149.0 ± 0.3
ALM040280	198.9 ± 0.4
ALM040287	251.1 ± 0.5
ALM040284	318.6 ± 0.6
ALM040278	441.7 ± 0.9



STATE OF MISSOURI

DEPARTMENT OF HEALTH AND SENIOR SERVICES BREATH ALCOHOL PROGRAM



PERMIT TYPE II **BRYAN J REID**

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

DATAMASTER, INTOX EC/IR II

for the determination of the alcoholic content of blood from a sample of expired air. Permit issued under the provisions of sections 577.020 through 577.041, RSMo and 306.111 through 306.119 RSMo.

DATE	white
	DIRECTOR OF STATE PUBLIC HEALTH LABORATORY
NUMBER 240316	Dal Vasterly
EXPIRES 7/30/2016	Zioca Vastacen
40 580-0771 (8-10)	DIRECTOR OF DEPARTMENT OF HEALTH AND SENIOR SERVICES

LAB-4 (R6-10)

